



CCLRC

Rutherford Appleton Laboratory

Laser status

CTF3 Collaboration meeting

CERN 23-25 November 2004

**Marta Divall, Graeme Hirst, Emma Springate, Ian Ross, Guy Suberlucq,
Roberto Losito**

M.Divall

CLF

Reminder to the PILOT test results, Pilot to CTF3



- Milestones and deliverables
- Development programme
- Stability
- Steady-state operation
- Oscillator, preamplifier
- Thermal management
- Amplifier design (gain calculations)
- High power pumping diodes
- Amplifier head design
- Active feedback stabilisation
- Summary

M.Divall

CLF

PILOT TEST

System layout



June/2002

Divall

CLF

PILOT to CTF3

Target parameters

Target Parameters	Units	PILOT (achieved)	CTF3
Charge / bunch	nC	0.07	2.33
Number of pulses	-	250	2332
Oscillator repetition frequency	GHz	0.25	1.5
Macro-pulse width	μ s	1	1.548
QE _{min}	%	4	3
Wavelength	nm	262	
E _{cathode} / pulse	nJ	8.75	370
Optical path transm.	%	52	70
IR/UV conversion eff.	%	1.5	7.4
Stabilization transm.	%		70
Repetition rate	Hz	5	5 (1-50)
E _{OUT} / pulse (Amplifier)	μ J	1.12	10
Pulse train mean power	kW	0.28	15
Average power (UV)	mW	1.4	116
E _{OUT} / pulse (oscillator)	nJ	0.4	6.7

M.Divall

CLF

Milestones

- Increase the power of existing oscillators (CCLRC-RAL, CNRS-LOA).
- Optimise the design of the amplifiers for required power output at minimum cost and complexity (CCLRC-RAL, CERN, CNRS-LOA).
- Efficiently convert the laser output wavelength to UV using new harmonic generator crystals (CCLRC-RAL, CERN, CNRS-LOA).
- Establish ultra-high stability including the use of charge and pointing stability feedbacks (CCLRC-RAL, CERN, INFN-LNF).

Deliverables:

- High power oscillator.
- Specific amplifiers.
- Specific frequency conversion stages.
- Test of feedback systems.

M.Divall

CLF

- Photo-cathode performance

- Oscillator

high cw power ~30W
1.5GHz mode-locked

- Amplifier

stable output pulse train
high efficiency (1 μ m IR laser)
good beam quality (lensing-effect)

- Feedback-control (<0.1%rms)

laser diode pump current (μ s)
fast optical gate (Pockels-cell)
accurate monitor (0.1%)

- Fourth harmonic generation

high efficiency
thermal management

- Optics

good beam quality
compensation for thermal lensing
high damage thresholds

- Diagnostics

stable operation
high resolution for stability measurements
1-50Hz qcw signal detection

CTF3 amplifiers

Stability

- Diode pumping
 - Steady-state operation
 - Feedback-system (Pockels cell, pumping current)
 - Saturated amplifiers
 - Saturated harmonic generation stages
-
- Temperature of the coolant (wavelength shift of diodes)
 - Pump current stability
 - Oscillator stability
 - Ringing in the Pockels cell(s)
 - Thermal effect in the harmonic generation stages

M.Divall


CLF

Pulse structure




1.5 GHz cw pulse train

X 300



Pumping pulse




333ps



Phase coding after
or before the amplifiers


X 5




Pumping pulse



1.5 GHz 3 GHz



Gating pulse





M.Divall

CLF

Steady-state operation

$$P_{out} - P_{in} = \eta_p P_{abs} - \frac{\pi D^2}{4} F_{sat} \ln G \frac{(1+B)}{\tau_{fl}}$$

Effect of seed modulation on the output






Fast feedback stabilisation is needed!!

M.Divall

CLF

Feedback stabilisation



0.25%rms stability over 150 μ s

M.Divall

CLF

~4 μ s response time
+/- 10% change can be compensated
~65 % efficient

Speed needs to be improved (faster electronics)
Built in compensation for nonlinearity
Tests on the real system

Oscillator and preamplifier

Parameter	unit	Specified in tender	HighQ
Laser gain medium		Nd:YLF	
Wavelength	nm	1047	
Type of operation		cw mode-locked	
Repetition rate	MHz	1499.28	
Average output power	W	> 0.2	>0.2 3W, 10W with preamplifiers
Pulse width (FWHM)	ps	< 10	
Timing jitter from an external 1.5 GHz RF source	ps	<±1	<±1 (typ.)
Polarization		Linearly polarized (1:500)	
Beam quality M^2		$TEM_{00} M^2 < 1.2$	$M^2 < 1.2$
Beam pointing stability	$\mu\text{rad}/^\circ\text{C}$	±25	±25
Beam size stability		5% rms size jitter	
Amplitude stability after a 1 h warm-up		< 0.2 % rms above the 100 kHz noise region <1% rms below the 100kHz noise region	
Warranty	year	1	1year/5000h of operation whichever is first
Repair and maintenance		Available for at least 5 years	

Tender has been sent out 28th of September 2004

Delivery will be February of 2005

The acceptance test will take place in CERN

The support will be valid in RAL after shipping

The company will provide an all optical timing jitter measurement

M.Divall

CLF

CTF3 amplifiers

High average power, thermal management

- Thermal lensing
- Fracture limit of YLF for rod geometry 22 W/cm

With **37 kW** peak pumping power the thermal load **~21 W/cm** assuming homogeneous pumping along the length of the rod.



Fracture limit measurements for etched rods
Redesign the system with less power in the final amplifier


Possible collaboration with Vienna University of Technology

- Development of etching technology for strengthening of LiYF₄ laser crystals
- Development of ion exchange technology for strengthening of LiYF₄ laser crystals
- Combination of etching and ion exchange techniques for strengthening of LiYF₄ laser crystals

M.Divall

CLF


New system layout



M.Divall

CLF

Calculated gain curves



M.Divall

CLF

Advantages

- Less power in the second amplifier; safer thermal conditions
- 50 Hz operation is possible
- Higher overall efficiency
- Smaller gain and higher saturation in the second amplifier
- Similar power level in the amplifiers allows easier thermal lens compensation
(The two crystals can be rotated by 90° to each other)

Disadvantages

- Two completely new amplifier
- Higher cost

New specification for the diodes

M.Divall

CLF

Pumping diodes

1st amplifier

15 kW peak power 400 μ s 50Hz pumping

Asking for

- 18 kW peak power for long lifetime 10^9 shots
- 2.5% duty cycle for flexibility in pumping time
- Homogeneous pumping along the length of the 7 cm rod
- 5 fold geometrical configuration

2nd amplifier

20 kW peak power 200 μ s 50Hz pumping

Asking for


- 22 kW peak power for long lifetime 10^9 shots
- 2% duty cycle for flexibility in pumping time
- Homogeneous pumping along the length of the 11 cm rod
- 5 fold geometrical configuration

Specification needs to be put down for diodes
Homogeneous pumping has to be investigated
Surveying on the market has to be carried out


M.Divall

CLF

The geometry




Astigmatism (As)



5kW
7.5 cm **Vertical stack**
Need bigger spacing between arrays

Vertical stack	Horizontal stack
Mixed polarization \perp and \parallel to c-axis	Only perpendicular Polarization to c-axis
Some parallel to c-axis \rightarrow higher absorption efficiency	More sensitive to Wavelength and so temperature
Too high power/length \rightarrow several stacks needed or bigger spacing	More flexible in power/length
The divergence angle is ideal for 1 cm diameter rod	Special microlens arrangement needed for 6-8° divergence
Complicated to calculate pumping distribution	Easy to calculate pumping distribution




5.5 kW
10.92 cm **Horizontal stack**
Lensed option down to 1°


M.Divall

CLF

Pumping distribution Measurements vs. calculations





Changing the position of the diodes




Measurements

Changing the temperature +/- 5°C




Fa2spec



Fa2spec

792nm perpendicular to c-axis
(arrays along the optical axis)
0.5cm diameter rod

MathCad
calculations



LasCad
calculations


M.Divall

CLF

Pumping diodes

- Lensed option is too expensive, and no guarantee for good performance
- To buy more arrays to fill the length of the rod is cheaper, than a custom bar spacing
- The estimated cost for the two amplifiers together
max. 400.000 Euro

The design calls for uniform pumping along the effective 7 cm of the rod, using either a single long stack or 2-3 shorter stacks. External lenses will assure efficient coupling of the pump into the rod.




The tender was out by 1st of November

M.Divall

CLF

The head design

- Simple water and electrical connections
- More accurate positioning of the focusing optics
- Easy assembly
- Possibility of rotating the rod in situ
- Alike design for the two amplifiers



M.Divall

CLF



Referees suggestions

- Utilizing the energy stored in the laser rods during the total pump-up phase in more efficient way would probably allow to reduce the required amount of diodes and therefore the total costs of the laser system. It would also solve the problems with fracture limit of the laser rods. This option should be taken into account because it is of large interest for other institutions.

With burst mode less pumping time is needed. However to get the energy per pulse out from the system (peak power) the same pumping power needed. The required duty cycle for the diodes is standard, and the lifetime will depend on the number of shots, and not on the energy/pulses. So cost-wise the qcw solution for the seed does not make a difference. However it would help with thermal management, but in the new design this problem is solved, and with thermal power/unit length being the same in both amplifiers the thermal lensing can be easily compensated. With steady-state operation in the case is slow input variations the output will remain the same, while in the case of burst mode operation this can cause very high overshoot (gain) at the front of the train. For safe operation the system requires sophisticated feedback stabilization and timing control.

- It may worthwhile to consider an oscillator running in burst-mode instead of CW mode.

This type of laser is not common on the market. Stable operation is guaranteed in a cw mode-locked laser.




IEEE JOURNAL OF QUANTUM ELECTRONICS, VOL. 34, NO. 10, OCTOBER 1998

Feedback-Stabilized Nd:YLF Amplifier System for Generation of Picosecond Pulse Trains of an Exactly Rectangular Envelope
Ingo Will, Armin Liero, Dieter Mertins, and Wolfgang Sandner

M.Divall

CLF

Summary



M.Divall

CLF